Medtronic

**ZEISS KINEVO®* Microscope Interface for
Stealth AXiS™ Surgical Systems with
Cranial Quick Reference Guide**

Contents

English	3
Dansk	16
Deutsch	29
Español	44
Suomi	59
Français	72
Hrvatski	85
Italiano	98
Nederlands	111
Norsk	125
Polski	138
Português	151
Português do Brasil	164
Русский	177
Svenska	191
Українська	203

EN

English

Description

The Medtronic Navigation microscope tracker is an accessory to the Stealth AXiS™ Surgical System. It consists of a tracker assembly that attaches to the ZEISS KINEVO microscope.

Intended use

The Medtronic Navigation microscope tracker is an accessory to the Stealth AXiS™ Surgical System. It is intended to be used with the ZEISS KINEVO microscopes and the Stealth AXiS™ Surgical System to enable microscope tracking positions during surgical procedures.

Indications for use

The Stealth AXiS™ Surgical System is intended for precise positioning of surgical instruments and as an aid for precisely locating anatomical structures in open, minimally invasive, and percutaneous procedures. The Stealth AXiS™ Surgical System is indicated for medical conditions in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure, such as the skull or vertebra, can be identified relative to a CT or MR based model, fluoroscopy images, or digitized landmarks of the anatomy.

Contraindications

None.

Material composition

Material contained in product that can cause an allergic reaction: Nickel

For additional materials of concern information such as REACH, CA Prop 65, or other product stewardship programs, go to www.medtronic.com/productstewardship.

Compatibility

This device is compatible with Medtronic Stealth AXiS™ computer-assisted surgery systems and Medtronic navigation instruments.

System level manuals

This device can be used with Medtronic computer-assisted surgery systems. See the Stealth AXiS™ system level manuals and clinical application for descriptions of patient groups, intended users, potential clinical benefits, side effects, and potential complications.

Qualified microscopes

When used with the Stealth AXiS™ Surgical System with Stealth AXiS™ Cranial clinical application, the Medtronic Navigation microscope tracker is qualified with the following ZEISS KINEVO microscopes: ZEISS KINEVO 900

Warnings

Warning: Inspect the microscope tracker assembly for damage prior to use. If visibly damaged or bent, do not use the microscope with the Stealth AXiS™ Surgical System. Do not attempt to repair it.

Warning: Use care when handling the microscope. Do not allow the microscope tracker to collide with anything.

Warning: If the microscope tracker assembly is damaged during the procedure, abort navigation and abandon use of the microscope with the Stealth AXiS™ Surgical System.

Warning: Frequently confirm localization accuracy and system responsiveness during navigation. If accuracy degrades, re-register the patient. Refer to the Stealth AXiS™ Cranial clinical application for instructions on how to re-register the patient.

Warning: The device is designed for re-use and is supplied non-sterile.

Warning: Use caution and exercise care in order to prevent injury during use. See the corresponding system level manuals and clinical application for further handling instructions.

Warning: To avoid potential exposure to blood-borne pathogens and chemicals, use appropriate Personal Protective Equipment when handling, processing, or disposing of Medtronic devices.

Warning: Confirm navigational accuracy via anatomical landmarks and accuracy checkpoints before navigation.

Precautions

Caution: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

Caution: Refer to the Stealth AXiS™ system level manuals and clinical application for additional instructions, warnings, and cautions.

Setting up equipment

1. Connect the Stealth AXiS™ Surgical System to a power source and turn it on.
2. Connect the ZEISS KINEVO microscope to a power source and turn it on.
3. Connect the cables.

Note: Do not connect to an external display device while performing a clinical procedure. When connecting or disconnecting an external monitor, all monitors on the system will go black for a few seconds while the video subsystem reconfigures.

Note: The ZEISS KINEVO microscope interface cable provides connectors to transmit video, data, and tracker information between the Stealth AXIS™ Surgical System and the ZEISS KINEVO microscope.

- a. If you are connecting the ZEISS KINEVO microscope and the Stealth AXIS™ Surgical System to the hospital network, use two Ethernet cables to make the following connections:

Note: The Ethernet cable transmits the microscope settings between the ZEISS KINEVO microscope and the Stealth AXIS™ Surgical System.

- Connect one Ethernet cable between the ZEISS KINEVO microscope and a wall jack.
 - Connect the other Ethernet cable between the Stealth AXIS™ Surgical System and a different wall jack.
- b. Attach the following connectors on the ZEISS KINEVO microscope interface cable to the ZEISS KINEVO microscope.
- First, insert the large, rectangular ODU™* connector **(1)** and turn the knob to secure.
 - Insert the video connector to one of the video in ports, Port 3 or Port 4 **(2)**.
 - If you are not using the hospital network, insert the network Ethernet connector to the network port **(1)**.
- c. Attach the following connectors on the ZEISS KINEVO microscope interface cable to the Stealth AXIS™ Surgical System.
- Connect the yellow latching push-pull LEMO™* connector to the yellow Microscope port on the instrument interface panel of the Stealth AXIS™ Surgical System **(3)**.
 - Connect the gray latching push-pull LEMO™* connector to the Instrument port 1 or 2 on the instrument interface panel of the Stealth AXIS™ Surgical System **(4)**.
 - If you are not using the hospital network, insert the network Ethernet connector to the network port.

Figure 1: Connect the cables

- 1 ZEISS KINEVO microscope tracker and network connections
- 2 ZEISS KINEVO microscope video connection
- 3 Microscope port (digital) on Stealth AXIS™ Surgical System
- 4 Instrument ports 1 and 2 on Stealth AXIS™ Surgical System

4. Drape the microscope using an appropriate ZEISS KINEVO microscope drape. Refer to the instructions for use that accompany the ZEISS KINEVO microscope for proper draping.
 - a. Pull the drape taut over the LED tracker on the microscope head.
 - b. Secure the drape in place with sterile rubber bands.

Note: Maintain visibility of the LEDs by making sure that they are clearly visible and unobstructed by the rubber bands and any wrinkles in the drape. Make sure that the LEDs are clearly visible and unobstructed by the rubber bands and any wrinkles in the drape.

5. Perform an autobalance of the ZEISS KINEVO microscope. Refer to the instructions for use that accompany the ZEISS KINEVO microscope.
6. Set up tracking.
 - a. In the **Equipment** task in the clinical application, make sure that the microscope is connected to the Stealth AXIS™ Surgical System.
If the microscope equipment card is not displayed, search for the microscope in the **Add Equipment** list and tap **[+]** next to the microscope equipment card.

If the microscope equipment card indicates that the microscope is not connected to the clinical application, refer to the “Troubleshooting” section.

Note: If you have two or more active instruments in a procedure, assign each instrument to an instrument port. See the clinical application instructions for use for information about how to assign an active instrument.

- b. Tap the Settings button and confirm that the IP address shown matches the IP address of the ZEISS KINEVO microscope. To obtain the IP address of the microscope, refer to the instructions for use that accompany the ZEISS KINEVO microscope.
- c. On the ZEISS KINEVO microscope, activate the Navigation (extended) feature.

Note: During the activation process, the ZEISS KINEVO microscope displays a “Connection Attempt” message where you select the correct Display Port. If you connected the video cable to Port 3, select **DisplayPort 3**; if you connected the video cable to Port 4, select **DisplayPort 4**. Then, under Navigation (extended), choose **Active** under the display port that you chose.

Note: After communication between the Stealth AXiS™ Surgical System and the ZEISS KINEVO microscope is established, a green bar appears on the equipment card. The Stealth AXiS™ Surgical System and the ZEISS KINEVO microscope each display the hash key (the unique identifier for the microscope) and the serial number of the other device.

- d. Position the camera so that both the microscope tracker and patient reference frame display in the tracking view on the clinical application. Make sure that the microscope tracker is the only tracked instrument in the camera’s field of view to avoid navigational interference.

Navigating

Warning: Before localizing, verify navigational accuracy of the system.

If navigational accuracy is affected, abort use of the system.

1. Verify navigational accuracy of the system by comparing points of interest on the clinical application to the corresponding points of the patient anatomy.
2. Localize the microscope.

To “localize” means to identify the location of the tip of a Stealth AXiS™ Surgical System probe on the preoperative patient images in the clinical application. When you use the microscope to localize, this “tip” is the anatomical point you observe at the intersection of the ocular crosshair in the focal plane of the microscope.

- a. Focus the microscope on the anatomy of interest.

- b. Verify the system responsiveness and accuracy by observing the updated localization site in the images on the clinical application. If localization seems inaccurate, follow the instructions in step 3.

Note: As you decrease magnification, the depths of field (DOF) increases. A larger DOF will result in objects further above and below the navigated focal point appearing in focus. To narrow the DOF around the navigated focal point, increase microscope magnification.

3. Focus on a reference frame divot.

Center the HUD-injected crosshair or converging lasers (if equipped) on a reference frame divot, increase magnification to the highest level, and focus the microscope on the bottom of the divot.

- Alternatively, you may focus on the center of the HUD-injected crosshair on the etched divot of the reference frame.
- Confirm that the banner message indicates that the distance to the divot is less than 2 mm.
- The localization site is at the intersection of the HUD-injected crosshair.

Warning: If inaccuracy is detected, the tracker assembly may have been moved or repositioned. Contact Medtronic Navigation to schedule a microscope install calibration. Calibration of the microscope tracker must be performed by a trained service person.


Heads-up display (HUD)

The Stealth AXIS™ Surgical System microscope interface can superimpose navigational information (like tractography, tumor boundaries, and surgical plans) on images that pass through one of the microscope oculars. To use HUD mode, first create a tract, model, or surgical plan and then use the microscope hand controls to select the desired HUD options.

Selecting HUD options using the microscope handgrip

- To view the selection menu through the ocular, press the Up, Down, Select, or Back button on the handgrip.
- To select the desired HUD option, highlight it in the menu and press the Select button.


Figure 2: Handgrip functions

- 1 Up
- 2 Back
- 3 Down
- 4 Select

HUD options

Figure 3: Microscope options

- **Plans & Annotations:** Injects the active surgical plan(s) and annotations, if they were added in the **Planning** or **Navigation** tasks on the clinical application.

Note: Only plans and annotations displayed in the views in the clinical application will be shown in the HUD.

- The plan's entry and target points display as circles and are joined by a solid line. The segment of the plan above the current focal plane of the microscope is semi-transparent. The segment of the plan below the current focal plane is opaque. The plan colors are similar to those displayed on the clinical application.
- The annotation displays as a sphere. An annotation above the current focal plane is semi-transparent. An annotation below the current focal plane is opaque and displays the annotation label. The annotation colors are similar to those displayed on the clinical application. If the annotation alert is turned on, the intersection of the alert and the focal plane is indicated using the same color as the annotation.

- **Models & Tracts:**

Note: Only models displayed in the 3D view in the clinical application will be shown in the HUD. Only tracts displayed in the views in the clinical application will be shown in the HUD.

- **Show Entire:** Shows the entirety of all models and tracts as viewed from the focal direction.
- **Show at Plane:** Shows the portion of the models and tracts as viewed from the focal direction. Only the intersection of the model and the plane is shown. You can adjust how many millimeters above and below the tracts will be displayed by updating the tract's card in the clinical application.
- **Show Distal to Plane:** Shows the portion of the models and tracts that are on the far side of the focal plane as if viewed from a focal direction.

- **Create Plan:** Creates a surgical plan in the Stealth AXiS™ Cranial clinical application based on the current position of the microscope. The plan is composed of the microscope's focal point (the plan target) and the microscope's trajectory (the plan trajectory).

Troubleshooting

Locate the issue in the following table and follow the action steps. If the symptom persists, contact technical services for assistance.

Issue	Action
The Equipment task shows no communication between the Stealth AXIS™ Surgical System and the microscope. Or The Equipment or Instruments task shows the microscope tracker as disconnected. Or The crosshairs and HUD menus do not display in the microscope ocular.	1. Make sure that the microscope is powered on. 2. Check the cabling path between the Stealth AXIS™ Surgical System and the microscope: a. Make sure the Stealth AXIS™ Surgical System and the ZEISS KINEVO microscope are both connected to the network or have a direct Ethernet connection. b. Make sure the ODU™ connector is attached to the microscope. c. Make sure that the yellow LEMO™ connector is connected to the yellow Microscope port on the instrument interface panel of the Stealth AXIS™ Surgical System. d. Make sure that the gray LEMO™ connector is connected to Instrument port 1 or 2 on the instrument interface panel of the Stealth AXIS™ Surgical System. e. If you are using multiple active instruments in a procedure, open the Instruments task and verify that the clinical application registers that the microscope is plugged in to the correct Instrument port. If the clinical application registers that the microscope is plugged into the wrong Instrument port, re-assign the microscope. 3. Go to the Equipment task on the Stealth AXIS™ Surgical System and verify that the equipment is set up correctly. a. Tap the equipment card for the ZEISS KINEVO microscope. The equipment card displays hash key information (the unique identifier for the microscope) and the serial number of the microscope. b. Follow the instructions on the equipment card to confirm the following settings on the ZEISS KINEVO microscope. ◦ Source is set to "Navigation" ◦ UDP Connection is "Enabled" ◦ The correct display port is selected ◦ Display Port connection is set to "Active"

Issue	Action
The microscope displays a red status bar in the Navigation task even though all cables are properly connected.	Position the camera so that both the microscope tracker and patient reference frame are in the field of view. For optimal positioning, align the camera to the forward right of the surgical area with the microscope tracker in the camera's field of view. Confirm that the drape is positioned properly over the tracker. If the problem persists, contact Medtronic technical services with the part number and serial number of the microscope tracker. The probe geometry error may be too high to use the microscope with the navigation clinical application.

Disconnecting the interface

Leave the microscope tracker permanently attached to the head of the microscope. This helps prevent invalidating the calibration.

Warning: Do not modify, disassemble, or remove the microscope tracker. Removing and reattaching the microscope tracker changes the position of the tracker in relation to the microscope. This invalidates the calibration of the tracker, and a new calibration is required to prevent inaccurate navigation. Calibration of the microscope tracker must be performed by a trained service person.

- 1. Disconnect the microscope communications cables from the Stealth AXiS™ Surgical System.
- 2. Disconnect the microscope communication cables from the microscope.
- 3. Coil and secure the cables on the microscope using Velcro™ ties.

Maintenance

- 1. Make sure that the microscope tracker is securely mounted to the microscope. If the microscope tracker is not securely mounted, contact Medtronic technical services.
- 2. Clean the microscope tracker.
 - a. Refer to the following table for compatible disinfectants and recommended chemistries.

Table 1: Compatible cleaner and disinfectant chemistries

Disinfectant type	Disinfectant level
Quaternary ammonium germicidal detergent solution	Low-level
Phenolic germicidal detergent solution	Low-level
Isopropyl alcohol	70%
Sodium hypochlorite	Low level 5.25% to 6.15%

- b. Wipe and clean the external surfaces of the microscope tracker with a moistened low-level disinfectant wipe per the disinfectant manufacturer's instructions. If a premoistened wipe is not available, use a clean wipe and moisten with the disinfectant. Do not over saturate the wipe so much that liquid can enter any openings.

Warning: Prevent fluid from entering any openings on the microscope tracker and cable.

- c. Wipe the external surfaces with a clean wipe moistened with water or 70% isopropyl alcohol to remove disinfectant residue on the microscope tracker.
- d. Allow the microscope tracker to air dry completely before reusing.
- e. Visually inspect the microscope tracker for cleanliness. If visible soil or residue remains, repeat the steps to clean the microscope tracker assembly.

3. Visually inspect the microscope tracker for damage.

Make sure that the microscope tracker does not have any of the following indications of damage:

- Illegible etchings or markings
- Visible cracking or bending
- Discoloration, corrosion, stains or rust

If visibly damaged, the microscope tracker is at the end of its useful life and should not be used with the Stealth AXiS™ Surgical System.

4. Clean the microscope interface cable.

- a. Refer to the following table for compatible disinfectants and recommended chemistries for the microscope interface cable.

Table 2: Microscope interface cable disinfectants

Disinfectant type	Disinfectant level
Quaternary ammonium germicidal detergent solution	Low-level
Phenolic germicidal detergent solution	Low-level
Isopropyl alcohol	70%
Sodium hypochlorite	Low level 5.25% to 6.15%

- b. Power down and unplug the system before performing cleaning or disinfection procedures.
- c. Unplug external cables.
- d. Clean the cable by wrapping a low-level disinfectant wipe completely around the cable starting at one connector.

Warning: Do not allow liquid to enter the interior surfaces of the connector.

- e. Wipe the cable from end to end by pulling the cable through the wipe. Ensure contact with the connector, bend relief, and cable.

- f. If disinfection is desired, make sure that the surface to be disinfected remains in contact with the disinfectant for the specified contact time as directed by the manufacturer.
 - The contact time refers to the minimum amount of time the surface needs to remain visibly wet.
 - Use additional wipes, if needed, to maintain continuous contact time with the disinfectant.
- g. Examine the surfaces for visible soil. If soil is present, repeat cleaning.
- h. Wipe the external surfaces with a clean wipe moistened with water or 70% isopropyl alcohol to remove disinfectant residue.

Specifications

Table 3: Microscope tracker specifications

Operating temperature	18°C to 30°C (64°F to 86°F)
Operating Relative humidity	30% to 75% Non-condensing
Storage temperature	10°C to 40°C (50°F to 104°F)
Storage Relative humidity	10% to 80% Non-condensing

Storage

Ensure equipment and packaging are dry before storing. Store in dry, clean conditions at ambient room temperature such that the equipment and packaging are not compromised.

Disposal

When the device shows visible signs of wear, it is at the end of its useful life. Dispose of it in accordance with your facility’s procedure for medical devices and national regulations.

Contact information

If a serious incident occurs in relation to the use of this device, report it to Medtronic technical services. If this incident occurs in the European Union, also report it to the competent authority in the Member State where the incident occurred.

Assistance

For questions about the device, or to order accessories, contact your Medtronic Navigation representative or call Medtronic Navigation at +1 800 595 9709 in the USA or +1 720 890 3160 for technical support.

Medtronic and Medtronic logo are trademarks of Medtronic. ®*/™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

This product may be protected by US Patents listed at: medtronic.com/patents.

Symbol definitions

The following symbols may appear on device packaging:

Rx Only

Caution: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

Do not dispose of this product in the unsorted municipal waste stream. Dispose of this product according to local regulations. See <http://recycling.medtronic.com> for instructions on proper disposal of this product.

China RoHS compliant: Environmental protection use period of 50 years.

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All Rights Reserved
©2026 Medtronic Navigation, Inc.
Medtronic Navigation, Inc.
Medtronic B.V.